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Dear Customer, Thank you for trusting our brand and buying Evorei BOOST 3.

Before using the device, please follow the user manual carefully.

01. SAFETY

- 1. Children of 8 years and older and people with reduced motor, mental or cognitive abilities can use the device only under guardian supervision.

02. ORAL HYGIENE IRRIGATOR Designed to clean the spaces between permanent braces, crowns, and dental bridges.

- 04. DEVICE DESCRIPTION 1. Tank lid with nozzle holders 2. Water tank (600ml)

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- 13. TROUBLESHOOTING PROBLEM Water leaking between the jet tip and the irrigator handle.

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Importer: KGK Trade sp. z o.o. sp. k. ul. Ujastek 5B 31-752 Kraków Polska

Wyrodkowano w P.R.C.

15. WARRANTY/ CLAIMS This product is covered by a 24-month warranty.

Complaints should be submitted with the complaint form which can be found at http://reklamacje.kgktrade.pl/

Contact details and the address of the service can be found on the website: www.evorei.pl

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste.

Importer: KGK Trade sp. z o.o. sp. k. Ujastek 5B 31-752 Cracow Poland

Made in P.R.C.